SC2Flu Triplex Fast Test (Colloidal Gold)  
For SARS-CoV-2 & influenza A/B antigens

*A POC test*  
*CIBG reg*

# SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD) FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS

## Specification

Component	Specification	Catalogue number			
		MU035A02	MU035A20	MU035A05	MU035A25
		Quantity	Quantity	Quantity	Quantity
<b>Test cartridge*</b>	1 per pouch	2	20	5	25
<b>Desiccant*</b>					
<b>Nasal swab</b>	1 per bag	2	20	5	25
<b>Sample extraction buffer</b>	400 µL/vial	2	20	5	25
<b>Instructions for use documents</b>	1 each - Instructions for use 1 each - Quick reference instruction card				

**\*Desiccant are packaged inside the test cassette pouch.**

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

## INTENDED USE

*The SC2Flu Triplex Fast Test (Colloidal Gold) For SARS-CoV-2 & influenza A/B antigens* is a colloidal gold-immunoassay intended to simultaneously qualitatively detect and differentiate the nucleocapsid protein antigen(s) of 2019 novel coronavirus (SARS-CoV-2), influenza A virus and influenza B virus from individuals suspected of respiratory viral infection in anterior nasal swab (ANS) specimen from individuals aged 14 years and older, or children aged 2-14 by their parents/guardians, or for prescription use only tests for symptomatic individuals who are suspected of COVID-19 or influenza disease by a healthcare provider, or individuals with or without symptoms or other epidemiological reasons to suspect these viral infection.

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

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
Professional use; POC test  
Home use

One Strip:  
Single sampling

C, control line; A, influenza A test line;  
B, influenza B test line; S, SARS-CoV-2 test line

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

## WHAT ARE THE ADVANTAGES?

**Fast! Result in 15 min or less**











POC test by Liberal nurse possible

Self test: Prescription or OTC

# Interpretation of the test results

Valid test		Invalid test
Image	Interpretation	Image
	Influenza A positive	
	Influenza B positive	
	SASR-CoV-2 positive	
	SASR-CoV-2/Influenza A pos	
	SASR-CoV-2/Influenza B positive	
	Influenza A/Influenza B positive	
	All negative	

**SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS  
CIBG/CE CERTIFIED**

**SAMPLE TYPE EQUIVALENCY**

Sample	Spiked virus	N	RSV Ag positive	SC2 Ag positive	FluA Ag positive	FluB Ag positive	Negative
NPS	500 TCID <sub>50</sub> /mL	10	10/10 (100%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	0/10 (0%)
OPS		10	10/10 (100%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	0/10 (0%)
Saliva		10	10/10 (100%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	0/10 (0%)
VTM		10	10/10 (100%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	0/10 (0%)
NPS	2500 TCID <sub>50</sub> /mL	11	11/11 (100%)	11/11 (100%)	11/11 (100%)	11/11 (100%)	0/11 (0%)
ANS		11	11/11 (100%)	11/11 (100%)	11/11 (100%)	11/11 (100%)	0/11 (0%)
Saliva		11	11/11 (100%)	11/11 (100%)	11/11 (100%)	11/11 (100%)	0/11 (0%)
VTM		11	11/11 (100%)	11/11 (100%)	11/11 (100%)	11/11 (100%)	0/11 (0%)
NPS	Saline	21	0/21 (0%)	0/21 (0%)	0/21 (0%)	0/21 (0%)	21/21 (100%)
ANS		21	0/21 (0%)	0/21 (0%)	0/21 (0%)	0/21 (0%)	21/21 (100%)
Saliva		21	0/21 (0%)	0/21 (0%)	0/21 (0%)	0/21 (0%)	21/21 (100%)
VTM		21	0/21 (0%)	0/21 (0%)	0/21 (0%)	0/21 (0%)	21/21 (100%)

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

## WHAT ARE THE ADVANTAGES?

Super sensitive!

Virus	Matrix	Limit of Detection
SASR-CoV-2	Anterior Nasal Swab	78 TCID <sub>50</sub> /mL
Influenza A	Anterior Nasal Swab	156 TCID <sub>50</sub> /mL
Influenza B	Anterior Nasal Swab	156 TCID <sub>50</sub> /mL

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

## High-dose Hook Effect

Virus	Matrix	High-dose Hook
SARS-COV-2	Anterior Nasal Swab	$4 \times 10^5$ TCID <sub>50</sub> /mL
Influenza B virus	Anterior Nasal Swab	$4 \times 10^5$ TCID <sub>50</sub> /mL
Influenza A virus	Anterior Nasal Swab	$4 \times 10^5$ TCID <sub>50</sub> /mL

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

Viral Ag tested swab specimens from healthy individuals

viral antigen tested		Nucleic Acid Tested		Total
		RNA pos	RNA neg	
Ag neg	SC2 Neg	0	151	151
	FluA Neg	0	143	143
	FluB Neg	0	162	162

**Negative percent agreement (NPA):**  
SC2 Ag testing NPA:  $151/151 = 100.00\%$  [97.52% ~ 100.00%];  
FluA Ag testing NPA:  $143/143 = 100.00\%$  [97.38% ~ 100.00%];  
FluB Ag testing NPA:  $162/162 = 100.00\%$  [97.68% ~ 100.00%].

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
CIBG/CE CERTIFIED

Sensitivity study with Nucleic Acid Testing

Viral antigen tested		Nucleic Acid Tested			Total
		SC2 RNA positive	FluA RNA positive	FluB RNA positive	
SC2 Antigen	Positive	32	0	0	32
	Negative	2	0	0	2
FluA Antigen	Positive	0	31	0	31
	Negative	0	1	0	1
FluB Antigen	Positive	0	0	39	39
	Negative	0	0	1	1
Total		34	32	40	106

**Positive Percent Agreement (PPA):**  
SC2 Ag testing PPA:  $32/34 = 94.12\%$  [80.32 ~ 99.28%];  
FluA Ag testing PPA:  $31/32 = 96.88\%$  [84.26 ~ 99.45%];  
FluB Ag testing PPA:  $39/40 = 97.50\%$  [87.11 ~ 99.56%].

**SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS  
CIBG/CE CERTIFIED**

# Clinical Trial Protocol in USA

## Clinical trial targets

*The sample matrix tested for both the candidate test and comparator test:*

**Anterior Nasal Swab (ANS)**

# Clinical Trial Protocol in USA

## Clinical trial targets

1. POC professional test;
2. Prescription home test;
3. OTC, home test;
4. Asymptomatic screen;
5. Confirmative test

## Clinical Trial Protocol: Test device determined

Item	Candidate Device	Comparator Device 2
Device Name	<i>SC2Flu Triplex Fast Test (Colloidal Gold) For SARS-CoV-2 &amp; Influenza A/B antigens</i> CIBG/CE CERTIFIED	cobas SARS-CoV-2 & Influenza A/B
Manufacturer	AMPER INC	Roche Molecular Systems, Inc.
Indication	to simultaneously qualitatively detect and differentiate the nucleocapsid protein antigen(s) of 2019 novel coronavirus (SARS-CoV-2), influenza A virus and influenza B virus from individuals suspected of respiratory viral infection in anterior nasal swab (ANS) specimen from individuals aged 14 years and older, or children aged 2-14 by their parents/guardians, or for prescription use only tests for symptomatic individuals who are suspected of COVID-19 or influenza disease by a healthcare provider, or individuals with or without symptoms or other epidemiological reasons to suspect these viral infection.	An automated multiplexed real-time RT-PCR assay intended for simultaneous qualitative detection and differentiation of SARSCoV-2, influenza A virus, and/or influenza B virus RNA in healthcare provider-collected nasal and nasopharyngeal swab specimens and self-collected nasal swab specimens (collected in a healthcare setting with instruction by a healthcare provider) from individuals suspected of respiratory viral infection consistent with COVID-19 by their healthcare provider. 1 Emergency use of this test is limited to authorized laboratories.
Authorized Laboratories	<ol style="list-style-type: none"> <li>1. POC professional test.</li> <li>2. Prescription home test.</li> <li>3. OTC, home test.</li> <li>4. Asymptomatic screen.</li> <li>5. Confirmation test.</li> </ol>	Laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet requirements to perform moderate or high complexity tests.
Date issue	Candidate	September 30, 2021
Kits needed	10,000 devices	10,000 SARS-CoV-2 Test Kits, 10,000 nucleic acid isolation & purification kits, 10,000 nasal swab collectors with VTM

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We, the manufacturer, here with declare that the product(s)

Product Name	<b>SC2Fu Triplex Fast Test (Colloidal Gold) For SARS-CoV-2 &amp; influenza A/B antigens</b>	Specification	<b>2 tests/box, 20 tests/box</b>
Intended Use	<i>The SC2Fu Triplex Fast Test (Colloidal Gold) For SARS-CoV-2 &amp; influenza A/B antigens is a colloidal gold-immunoassay intended to simultaneously qualitatively detect and differentiate the nucleocapsid protein antigen(s) of 2019 novel coronavirus (SARS-CoV-2), influenza A virus and influenza B virus from individuals suspected of respiratory viral infection in anterior nasal swab (ANS) specimen from individuals aged 18 years and older, or children aged 6-18 by their parents/guardians, or for prescription use only tests for symptomatic individuals who are suspected of COVID-19 or influenza disease by a healthcare provider, or individuals with or without symptoms or other epidemiological reasons to suspect these viral infection.</i>		
Classification	Others		

Conformity Assessment Route: **MDD98/79/EC Annex III.**

Applicable Standards:

<i>ISO 13485:2016</i>	<i>EN ISO 18113-2:2011</i>	<i>EN 13612:2002</i>
<i>ISO 14971:2019</i>	<i>EN 13641:2002</i>	<i>ISO 23540:2015</i>
<i>EN ISO 18113-1:2011</i>	<i>ISO 15229-1:2016</i>	<i>EN 62566-1:2015</i>

**CE**

We, the manufacturer, here declare with sole responsibility that our product(s) mentioned above meet(s) the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices. We agree to develop, implement and maintain a documented post-production monitoring process.

	
CEO, Dr. HUI-YU LIU, MD, PhD	Date & Place: July 1, 2021 @ NC, USA

Amper

## COMPANY PROFILE

- The AMPER is a high-tech company dedicated to the development, production, and sales of rapid IVD test devices. The company obtained the manufacture licenses from NMPA in 2010 and gained the "ISO13485:2016" certification. It is a "National High-tech Enterprise" and a "Provincial High-tech Enterprise". It has accumulated plenty of experience in colloidal gold technology from R&D to production. Our colloidal gold technology platform covers many products.

## OEM Manufacturer

GMP workshop: **2500 m<sup>2</sup>**

Assembly line: **9 lines**

employees: **200**

Capacity: **300,000 devices daily.**

CE certificates: **3 SARS-CoV-2**  
related Products with CE mark

Fully automatic manufacture line

*SC2FLU TRIPLEX FAST TEST (COLLOIDAL GOLD)  
FOR SARS-COV-2 & INFLUENZA A/B ANTIGENS*  
**CIBG/CE CERTIFIED**

- Manufacturer:

**AMPER, INC**

2500 Gateway Centre Blvd, Suite 400

Morrisville, NC 27560, USA

- email: [info@amperbio.com](mailto:info@amperbio.com)
- Web: [www.amperbio.com](http://www.amperbio.com)