


SOLMIRA®

CE 2934 IVD

SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Self-Test)



4 IN 1

Clinical signs and symptoms of respiratory viral infection due to SARS-CoV-2, influenza A/B and RSV can be similar. It is hard to make an accurate diagnosis based on symptoms alone.

- Differentiates between respiratory diseases
- Reduce the risk of a misdiagnosis
- Helps identify cases of co-infection

Cost-Effective

1 sample, 4 test results, 15 minutes to diagnosis

CE Mark for Self-Testing

CE self-testing certificated, Belgium AFMPS & Spain AEMPS Listed

Comprehensive

Effective for mutants Alpha, Beta, Gamma, Delta and Omicron of SARS-CoV-2; strains h3n2 and h1n1 of influenza A; strains victoria and yamagata of influenza B; types A and B of RSV

Self-Testing

Suitable for non-medical staff, no need for extra equipment.

Room Temperature

Transport, store and operate at room temperature

High Performance

	Sensitivity	Specificity
SARS-CoV-2	92.93%	100.00%
Influenza A	92.04%	100.00%
Influenza B	90.91%	100.00%
RSV	95.45%	100.00%


Components of 1 test




Characteristics

Method	Colloidal Gold Chromatographic Immunoassay
Reaction Time	15 mins
Shelf Life	18 months
Storage Temperature	2-30°C (Room temperature)
Sample Types	Nasal swab
Specification	1 Test/box, 2 Tests/box, 5 Tests/box


Brief Operation Procedure




1 Twist off the cap of the diluent tube and remove the inner blue stopper. Insert the diluent tube into the hole of the kit or use other items to hold the diluent tube in place.



2 Insert the nasal swab into the left nostril to a depth of 2.5 cm (1 inch) from the edge of the nostril. Rotate the nasal swab on the nostril wall (mucous membrane) 5 times to ensure adequate sampling. Repeat the process in the right nostril with the same nasal swab.



3 Place the swab sample into the tube, then break the swab at the swab node and leave the lower half in the tube. Close the cap. Squeeze the swab 10 times in the test tube.



4 Unscrew the terminal at the top of the cap. Each sample well of the test card requires 2 drops (about 60 µL) of the treated sample solution. Keep the test card at room temperature for 15 minutes to observe the test results.



Shenzhen Microprofit Biotech Co.,Ltd.

Tel: +86-755-61688835

Fax: +86-755-61688111

E-mail: bio@microprofit-bio.com


Website: www.microprofit-bio.com



Shenzhen Microprofit Biotech Co.,Ltd.

Add: Rm 405,406 Zone B/4F, Rm.205,206-1,207, West Side of Zone B/2F, Haowei Building, NO.8 Langshan 2nd Road, Songpingshan, Songpingshan Community, Xili Street, Nanshan District, Shenzhen, P.R. China

CE Certificate



Belgium AFMPS & Spain AEMPS Listed

Belgium

Lijst van aanbevolen SARS-CoV-2 antigen zelftesten
Liste des autotests antigéniques du SARS-CoV-2 recommandés
List of recommended SARS-CoV-2 antigen self-tests
version/version 29.07.2022

Fabrikant/Fabricant/Manufacturer	gemachtigde/representant autorise/authorised representative	Naam van de test/Nom de test/Name of the test	Referentie/Référence/Reference (marques commerciales/handelsmerken/trademarks)
Shenzhen Microprofit Biotech Co., Ltd. (CN)	CMC MEDICAL DEVICES & DRUGS, S.L. (ES)	SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay)	MF-71 (MF71-1; MF71-2; MF71-5) (Fluorecare)

*Les versions de ces documents sont les plus récentes reçues par l'AFMPS (des mises à jours ont pu être réalisées)
 Deze versies van de documenten zijn de meest recente die het FAGG ontvangen heeft (het zou kunnen dat er updates zijn geweest)
 The versions of these documents are the most recent received by the FAMHP (updates may have been made)*

Spain

Fluorecare Kit de Prueba Combinada para el Antígeno del SARS-CoV-2, la Gripe A/B y el VSR	Shenzhen Microprofit Biotech Co., Ltd. (China)	MF-71-1 (1 test); MF-71-2 (2 test); MF-71-5 (5 test)	Nº 2934	NASAL
---	--	--	---------	-------

About Us

Founded in 2009, Microprofit Biotech is a national high-tech enterprise from China, which is specialized in medical devices and reagents in the field of in-vitro diagnosis (IVD) with the idea of POCT (point-of-care-testing). fluorecare®, the official brand of Microprofit Biotech, dedicated to entire products of point-of-care (POC) tests based on immunofluorescence quantitative assay platform.



With professional R&D, manufacture, marketing, sales and service team, Microprofit Biotech sold products to nearly 100 countries&areas and has won a good reputation around the world. Microprofit Biotech strictly complies with national and international regulations and standards and got many certification, such as CE certificates, FDA certificates, FSC certificates and ISO 13485:2016 quality system certification.

